



## PRESTAN AED ULTRA TRAINER

**SKU:** PP-AEDUT-101

The AED for training purposes. PRESTAN is the First and Only Universal AED Trainer with **2015 Guidelines Compliant and Upgradable**.

The PRESTAN AED UltraTrainer is equipped with a dual graphics pad and Prestan's unique Pad Sensing System, which automatically detects when each pad is placed on the manikin. This AED trainer can accurately simulate a live AED, allowing for a more realistic training class.

The pads are durable, lasting up to 100 applications, and do not leave behind that annoying sticky residue. Replacement is easy and affordable.

The PRESTAN AED UltraTrainer is the first and only universal trainer to have a child button that best simulates modern live AEDs, giving your students a realistic training experience.

The easy to use audio menu allows you to customize your training with options of automatic or semi-automatic shock, 5 training scenarios, and the ability to turn on or off the compression metronome and ventilation command.

Features:

- PRESTAN is the First and Only Universal AED Trainer with **2015 Guidelines Compliant and Upgradable**
- **Play/Pause Button**
- **5 Training Scenarios**
- Dual Graphic (Adult/Child) Training Pads
- Pad Sensing System
- Preconnected Pads
- Child Training Button
- Child AED specific voice prompts
- Automatic or Semi-Automatic shock simulation
- CPR Coaching, "Give Breaths" voice prompt
- Turn on or off compression metronome
- Turn on or off "Give Breaths" prompts
- Low Battery warning

---

Download PDF:



[Read More](#)



## PRESTAN PROFESSIONAL AED TRAINER PLUS

SKU: PP-AEDT2-101

### QUALITY / DURABILITY / MOST SIMILAR TO A LIVE AED

- Leading-edge AED Trainer that most closely resembles brand AED Trainers in size and weight for a realistic training experience.
- It features customizable options to match training environments and includes a child button.
- Clear confident voice prompts. Also, CPR prompts to help the student keep pace with chest compressions.
- Easy to replace modules containing current CPR guidelines and 2 languages.
- A remote is optional and available for instructors who prefer this feature. A single remote will control several AED Trainers within its range.



Change the Language with a simple press of a button



Scenario button can quickly change the shock sequence (no need for manual)



Three volume options: loud, medium, soft



Play and Pause button for those teachable moments

#### Features:

- **Compliant with current 2015 AHA Guidelines**
- Fast and Easy as well as Update
- Adult/Child selection buttons
- Training pads with dual adult/child graphics
- Brightly illuminated shock button that flashes when the shock is advised
- 5 preprogrammed training scenarios with Clear and calm voice prompts
- Chest compression rate metronome with audible count reminders
- Automatic or semi-automatic shock delivery
- Programming option to either include or not include time for breaths after the compression cycle
- Available remote control (**sold separately**)
- Long battery life with automatic shut-off
- Replaceable language module
- Sturdy, flip-up lid

Download PDF:



[Read More](#)



## PRESTAN PROFESSIONAL AED TRAINER REMOTE

**SKU:** PP-AEDT-100-R

To use with the Prestan Professional AED trainer. Remote for Prestan Professional AED Trainer by using 2 AA batteries. Remote control lets the Instructor select a scenario, pause to give instruction, and then resume the scenario in progress.

Features:

- Set-up buttons
- Unit control buttons:
- Simulate loose pad
- Initiate "Press Deeper" during CPR
- Pause
- Play
- Send scenario

[Read More](#)



## PRESTAN PROFESSIONAL AED TRAINER PLUS PADS

**SKU:** PP-APAD2-4

PRESTAN Professional's high-performance training pads feature clear anatomical placement graphics and a long-life adhesive.

The pads are made of the same foam material used in live AED pads and conform well to all CPR manikins on the market today. The pad's silicone-based adhesive, coupled with their foam structure, ensure easy application and removal from your manikin.


The manufacturer states that they are good for at least 25 sticks and release AED training skill sessions, but we've used them for many more than that! The AED training pads ship in a zippered plastic bag, but once connected, stay clean and protected in the unit's training pad case.


**\*Pads are for training purposes only and cannot be used to deliver a defibrillation shock.**

Features:

- Fits Prestan AED Trainers
- Includes pad sensor
- Easy and clean removal from manikins
- Works with any brand manikin
- Latex-free silicone adhesive

[Read More](#)





## PRESTAN AED ULTRATRAINER PADS (4 PAIRS IN A PACK)

**SKU:** PP-UTPAD-4

Prestan AED UltraTrainer™ Adult/Child Replacement Training Pad Set, 4-Pack (8 pads total).

Prestan AED UltraTrainer pads are durable and long-lasting and feature an exclusive pad sensing system that allows the trainer to determine the pads are attached to the manikin. The silicone-based adhesive provides easy application and removal, helping prolong the life of the manikin skin.

Features:

- Dual Graphic (Adult/Child) Training Pads
- Pad Sensing System
- Preconnected Pads

[Read More](#)



