



PRODUCT CATALOG

generated from prestan.asia



TABLE OF CONTENTS

Parts and Accessories 2

PRESTAN CPR Manikins 12

PRESTAN AED Trainers 60



PRESTAN PROFESSIONAL ULTRALITE 50PCS FACE-SHIELD / LUNG BAG

SKU: PP-ULB-50

The face shield lung bag delivers a realistic experience of a full chest rise during CPR training. These lung bags fit the Prestan Ultralite series CPR manikins only.

This pack includes 50 lung bags. Insertion of the one-piece lung bag design is a breeze! Fold it up, slide it in and you are ready for training.

Features:

- Easy to insert
- Work with Prestan Ultralite CPR Manikins
- No tools necessary to insert

[Read More](#)



PRESTAN PROFESSIONAL INFANT 50PCS FACE-SHIELD / LUNG BAG

SKU: PP-ILB-50

Replacement for Manikin Face Shield Lung Bags

Features:

- Prestan Manikin Replacement Face Shield Lung Bags
- Clear plastic inflatable lungs
- Easy to insert
- Disposable

[Read More](#)





PRESTAN PROFESSIONAL CHILD 50PCS FACE-SHIELD / LUNG BAG

SKU: PP-CLB-50

Replacement for Manikin Face Shield Lung Bags.

Features:

- Prestan Manikin Replacement Face Shield Lung Bags
- Clear plastic inflatable lungs
- Easy to insert
- Disposable

[Read More](#)



PRESTAN PROFESSIONAL CHILD 50PCS FACE-SHIELD

SKU: PP-CFS-50

In classes where there is more than one student per manikin, there is a need for a balance between protecting the student and giving the student a realistic experience.

This is a 50 pack of Filtered Face Shields for Prestan Professional CPR Manikins. If manikins are shared during your training class, you'll want to use a filtered face shield to minimize the risk of cross-contamination between students. The Prestan face shield makes it easy and since it's secured to the manikin, it won't slip. No need for adhesives! Your students will spend less time distracted by an unruly barrier and more time practicing their skills.

Every CPR/AED instructor wants the best outcome when and if the student has a real-life experience giving CPR. This face shield, placed over the manikin's face, allows each student to practice CPR, be protected, and have the benefit of a realistic experience of a full chest rise.

50 pack face shields for Prestan Child Manikins.

Features:

- Professional Child Face Shields by Prestan Products (50 Pack)
- Not made with natural rubber latex

[Read More](#)





PRESTAN PROFESSIONAL ADULT 50PCS FACE-SHIELD / LUNG BAG

SKU: PP-ALB-50

Manikin Replacement Lung Bags

Features:

- Prestan Manikin Replacement Lung Bags.
- Clear plastic inflatable lungs.
- Easy to insert.
- Disposable.

[Read More](#)



PRESTAN FEMALE ACCESSORY


SKU: PP-ULB-50-1

The PRESTAN Female Accessory™ (*patent-pending*) has been designed to build confidence in CPR training, improve CPR survival rates in women, encourage dialogue & breakdown barriers, and eliminate gender bias in CPR administration.

Securely fits both the PRESTAN Professional Adult/Series 2000 & Ultralite Manikins. Pair with the PRESTAN CPR Training Shirts for a full real-life training experience.

Available in packs of 4 and in medium or dark skin tones.

[Read More](#)





PRESTAN MASKS ADAPTERS

SKU: 10076-PPA

Use this affordable, disposable adapter in place of expensive mask valves for cost-effective training. It features an integral air filter that obstructs the flow of contaminants. Fits all standard rescue masks (15mm I.D. or 22mm O.D.).

***For training purposes only.**

Features:

- It features an integral air filter that obstructs the flow of contaminants.
- Fits all standard rescue masks
- Rescue mask not included

[Read More](#)



PRESTAN PROFESSIONAL AED TRAINER REMOTE

SKU: PP-AEDT-100-R


To use with the Prestan Professional AED trainer. Remote for Prestan Professional AED Trainer by using 2 AA batteries. Remote control lets the Instructor select a scenario, pause to give instruction, and then resume the scenario in progress.

Features:

- Set-up buttons
- Unit control buttons:
- Simulate loose pad
- Initiate "Press Deeper" during CPR
- Pause
- Play
- Send scenario

[Read More](#)





PRESTAN AED ULTRATRAINER PADS (4 PAIRS IN A PACK)

SKU: PP-UTPAD-4

Prestan AED UltraTrainer™ Adult/Child Replacement Training Pad Set, 4-Pack (8 pads total).

Prestan AED UltraTrainer pads are durable and long-lasting and feature an exclusive pad sensing system that allows the trainer to determine the pads are attached to the manikin. The silicone-based adhesive provides easy application and removal, helping prolong the life of the manikin skin.

Features:

- Dual Graphic (Adult/Child) Training Pads
- Pad Sensing System
- Preconnected Pads

[Read More](#)



PRESTAN PROFESSIONAL AED TRAINER PLUS PADS

SKU: PP-APAD2-4

PRESTAN Professional's high-performance training pads feature clear anatomical placement graphics and a long-life adhesive.

The pads are made of the same foam material used in live AED pads and conform well to all CPR manikins on the market today. The pad's silicone-based adhesive, coupled with their foam structure, ensure easy application and removal from your manikin.


The manufacturer states that they are good for at least 25 sticks and release AED training skill sessions, but we've used them for many more than that! The AED training pads ship in a zippered plastic bag, but once connected, stay clean and protected in the unit's training pad case.

***Pads are for training purposes only and cannot be used to deliver a defibrillation shock.**

Features:

- Fits Prestan AED Trainers
- Includes pad sensor
- Easy and clean removal from manikins
- Works with any brand manikin
- Latex-free silicone adhesive

[Read More](#)





PRESTAN FEMALE SERIES 2000 MANIKIN

SKU: PP-AFM-2000-1-MS

The same high quality, realistic Professional Female Manikin is now available with Advanced Feedback including Rate, Depth, Recoil, Ventilation, Hands-Off Time, and Chest Compression Fraction via a Bluetooth® enabled PRESTAN CPR Feedback App available on both Apple and Android devices.

PRESTAN PROFESSIONAL FEMALE SERIES 2000 MANIKIN

The same high quality, realistic Professional Female Manikin is now available with Advanced Feedback including Rate, Depth, Recoil, Ventilation, Hands-Off Time, and Chest Compression Fraction via a Bluetooth® enabled PRESTAN CPR Feedback App available on both Apple and Android devices.

Fully compliant with current industry guidelines, including the Heart and Stroke Foundation of Canada's February 2021 Use of Feedback Devices in Resuscitation Training Courses Policy and the most recent AHA Integrated Feedback Directive.

Available in medium or dark skin tones and in Single or 4-Pack configurations.

All PRESTAN CPR Training Manikins are Made in the USA.

<https://youtu.be/7BD8vKb-CLc>

ADDITIONAL DETAILS:





****REDESIGNED** PRESTAN CPR Feedback App**

Provides high quality, real time, advanced CPR feedback via a Bluetooth® connection with the ****NEW & IMPROVED**** PRESTAN Professional Adult Series 2000 Manikin.


Advanced CPR feedback (includes rate, depth, recoil, ventilation, hands-off time, and ***chest compression fraction***) is monitored, easily displayed and available for reporting allowing students to self-monitor performance and gain confidence during CPR training.

The app allows instructors to provide comprehensive CPR training while monitoring CPR performance for up to six (6) students at a time and is equipped with a reporting feature that allows for CPR training results to be downloaded directly from the app.

The ****REDESIGNED**** PRESTAN CPR Feedback App provides an

- ***improved user experience***
- ***easier setup*** - with no need to calibrate breaths
- ***more intuitive CPR feedback*** - with the addition of chest compression fraction
- ***expanded list of available languages***

Languages currently supported by the PRESTAN CPR Feedback App include German, English (Intl), English (US), Spanish, French, Italian, Dutch, Russian, Korean, Arabic, Malay, Thai, Vietnamese, Belarusian, Croatian, Czech, Danish, Finnish, Greek, Hungarian, Norwegian, Polish, Portuguese, Romanian, Slovak, Slovene, Turkish, Ukrainian, and Chinese.



Download the ***REDESIGNED*** PRESTAN CPR Feedback App by Scanning the QR Code Below



The PRESTAN CPR Feedback App is for Instrumented Directive Feedback Devices (IDFD) that meet the most recent requirements of the American Heart Association and the Heart and Stroke Foundation of Canada's February 2021 *Use of Feedback Devices in Resuscitation Training Courses Policy*.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. FCC ID# QOQ-BGM220S2 / FCC ID# 2ASY3AK2G-PA

***NEW* PRESTAN Professional Adult Series 2000 Manikin CPR Rate Monitor LED Sequence**

	Connected Mode	Disconnected Mode
Sleep Mode/No Power		
Smart Device Connected Mode (Cycles in sequence left to right)		n/a
1 st Full Compression, Not Connected Mode (Wake mode/LED Test Mode)	n/a	
CPM* < 60		
CPM* 60 – 79		
CPM* 80 – 99		
CPM* 100 - 120		
CPM* > 120		
Compressions Stopped , Not Connected Mode		
Compressions Stopped , Not Connected Mode Low Battery Mode		
Compressions Stopped , Connected Mode		

*CPM = Compressions per minute
 = Flashing

In combination with PRESTAN’s traditional CPR feedback light system installed in the shoulder, the Series 2000 manikin offers a unique, dual-purpose training solution providing instructors with an option to train with either advanced or basic feedback.

The visual feedback gives students a real-life feel of delivering 100-120 compressions per minute. Lights continually signal rate progression, confirming compressions are within the 100- 120 cpm range recommended by the Industry

Guidelines. The yellow light will blink if the 120 cpm upper limit is exceeded, warning students to slow down. Instructors can monitor several students quickly and easily.

In addition to the visual feedback, students will also hear an audible beep as the chest is pushed to the appropriate depth of 2.0 -2.4 inches allowing them to experience the true force needed to deliver real life chest compressions.

[Read More](#)



PRESTAN INFANT SERIES 2000 MANIKIN

SKU: PP-IM-2000-1-MS

The same high quality, realistic Professional Infant Manikin is now available with Advanced Feedback including Rate, Depth, Recoil, Ventilation, Hands-Off Time, and Chest Compression Fraction via a Bluetooth® enabled PRESTAN CPR Feedback App available on both Apple and Android devices.

PRESTAN PROFESSIONAL INFANT SERIES 2000 MANIKIN

The same high quality, realistic Professional Infant Manikin is now available with Advanced Feedback including Rate, Depth, Recoil, Ventilation, Hands-Off Time, and Chest Compression Fraction via a Bluetooth® enabled PRESTAN CPR Feedback App available on both Apple and Android devices.

Fully compliant with current industry guidelines, including the Heart and Stroke Foundation of Canada's February 2021 *Use of Feedback Devices in Resuscitation Training Courses Policy* and the most recent AHA Integrated Feedback Directive.

Available in medium or dark skin tones and in Single or 4-Pack configurations.

All PRESTAN CPR Training Manikins are Made in the USA.

<https://youtu.be/7BD8vKb-CLc>



ADDITIONAL DETAILS:

***REDESIGNED* PRESTAN CPR Feedback App**

Provides high quality, real time, advanced CPR feedback via a Bluetooth® connection with the ***NEW & IMPROVED*** PRESTAN Professional Adult Series 2000 Manikin.



Advanced CPR feedback (includes rate, depth, recoil, ventilation, hands-off time, and **chest compression fraction**) is monitored, easily displayed and available for reporting allowing students to self-monitor performance and gain confidence during CPR training.

The app allows instructors to provide comprehensive CPR training while monitoring CPR performance for up to six (6) students at a time and is equipped with a reporting feature that allows for CPR training results to be downloaded directly from the app.

The ***REDESIGNED*** PRESTAN CPR Feedback App provides an

- **improved user experience**
- **easier setup** - with no need to calibrate breaths
- more **intuitive CPR feedback** - with the addition of chest compression fraction
- **expanded list of available languages**

Languages currently supported by the PRESTAN CPR Feedback App include German, English (Intl), English (US), Spanish, French, Italian, Dutch, Russian, Korean, Arabic, Malay, Thai, Vietnamese, Belarusian, Croatian, Czech, Danish, Finnish, Greek, Hungarian, Norwegian, Polish, Portuguese, Romanian, Slovak, Slovene, Turkish, Ukrainian, and Chinese.



Download the ***REDESIGNED*** PRESTAN CPR Feedback App by Scanning the QR Code Below



The PRESTAN CPR Feedback App is for Instrumented Directive Feedback Devices (IDFD) that meet the most recent requirements of the American Heart Association and the Heart and Stroke Foundation of Canada's February 2021 *Use of Feedback Devices in Resuscitation Training Courses Policy*.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. FCC ID# QOQ-BGM220S2 / FCC ID# 2ASY3AK2G-PA

***NEW* PRESTAN Professional Adult Series 2000 Manikin CPR Rate Monitor LED Sequence**

	Connected Mode	Disconnected Mode
Sleep Mode/No Power		
Smart Device Connected Mode (Cycles in sequence left to right)		n/a
1 st Full Compression, Not Connected Mode (Wake mode/LED Test Mode)	n/a	
CPM* < 60		
CPM* 60 – 79		
CPM* 80 – 99		
CPM* 100 - 120		
CPM* > 120		
Compressions Stopped , Not Connected Mode		
Compressions Stopped , Not Connected Mode Low Battery Mode		
Compressions Stopped, Connected Mode		

*CPM = Compressions per minute
 = Flashing

In combination with PRESTAN’s traditional CPR feedback light system installed in the diaper, the Series 2000 manikin offers a unique, dual-purpose training solution providing instructors with an option to train with either advanced or basic feedback.

The visual feedback gives students a real-life feel of delivering 100-120 compressions per minute. Lights continually signal rate progression, confirming compressions are within the 100- 120 cpm range recommended by the Industry

Guidelines. The yellow light will blink if the 120 cpm upper limit is exceeded, warning students to slow down. Instructors can monitor several students quickly and easily.

In addition to the visual feedback, students will also hear an audible beep as the chest is pushed to the appropriate depth of 1.5 inches allowing them to experience the true force needed to deliver real life chest compressions.

[Read More](#)



PRESTAN CHILD SERIES 2000 MANIKIN

SKU: PP-CM-2000-1-MS

The same high quality, realistic Professional Child Manikin is now available with Advanced Feedback including Rate, Depth, Recoil, Ventilation, Hands-Off Time, and Chest Compression Fraction via a Bluetooth® enabled PRESTAN CPR Feedback App available on both Apple and Android devices.

PRESTAN PROFESSIONAL CHILD SERIES 2000 MANIKIN

The same high quality, realistic Professional Child Manikin is now available with Advanced Feedback including Rate, Depth, Recoil, Ventilation, Hands-Off Time, and Chest Compression Fraction via a Bluetooth® enabled PRESTAN CPR Feedback App available on both Apple and Android devices.

Fully compliant with current industry guidelines, including the Heart and Stroke Foundation of Canada's February 2021 Use of Feedback Devices in Resuscitation Training Courses Policy and the most recent AHA Integrated Feedback Directive.

Available in medium or dark skin tones and in Single or 4-Pack configurations.

All PRESTAN CPR Training Manikins are Made in the USA.

<https://youtu.be/7BD8vKb-CLc>

ADDITIONAL DETAILS:

***REDESIGNED* PRESTAN CPR Feedback App**

Provides high quality, real time, advanced CPR feedback via a Bluetooth® connection with the ***NEW & IMPROVED*** PRESTAN Professional Adult Series 2000 Manikin.


Advanced CPR feedback (includes rate, depth, recoil, ventilation, hands-off time, and **chest compression fraction**) is monitored, easily displayed and available for reporting allowing students to self-monitor performance and gain confidence during CPR training.

The app allows instructors to provide comprehensive CPR training while monitoring CPR performance for up to six (6) students at a time and is equipped with a reporting feature that allows for CPR training results to be downloaded directly from the app.

The ***REDESIGNED*** PRESTAN CPR Feedback App provides an

- **improved user experience**
- **easier setup** - with no need to calibrate breaths
- more **intuitive CPR feedback** - with the addition of chest compression fraction
- **expanded list of available languages**

Languages currently supported by the PRESTAN CPR Feedback App include German, English (Intl), English (US), Spanish, French, Italian, Dutch, Russian, Korean, Arabic, Malay, Thai, Vietnamese, Belarusian, Croatian, Czech, Danish, Finnish, Greek, Hungarian, Norwegian, Polish, Portuguese, Romanian, Slovak, Slovene, Turkish, Ukrainian, and Chinese.



Download the ***REDESIGNED*** PRESTAN CPR Feedback App by Scanning the QR Code Below



The PRESTAN CPR Feedback App is for Instrumented Directive Feedback Devices (IDFD) that meet the most recent requirements of the American Heart Association and the Heart and Stroke Foundation of Canada's February 2021 *Use of Feedback Devices in Resuscitation Training Courses Policy*.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. FCC ID# QOQ-BGM220S2 / FCC ID# 2ASY3AK2G-PA

***NEW* PRESTAN Professional Adult Series 2000 Manikin CPR Rate Monitor LED Sequence**

	Connected Mode	Disconnected Mode
Sleep Mode/No Power		
Smart Device Connected Mode (Cycles in sequence left to right)		n/a
1 st Full Compression, Not Connected Mode (Wake mode/LED Test Mode)	n/a	
CPM* < 60		
CPM* 60 – 79		
CPM* 80 – 99		
CPM* 100 - 120		
CPM* > 120		
Compressions Stopped , Not Connected Mode		
Compressions Stopped , Not Connected Mode Low Battery Mode		
Compressions Stopped, Connected Mode		

*CPM = Compressions per minute
 = Flashing

In combination with PRESTAN’s traditional CPR feedback light system installed in the shoulder, the Series 2000 manikin offers a unique, dual-purpose training solution providing instructors with an option to train with either advanced or basic feedback.

The visual feedback gives students a real-life feel of delivering 100-120 compressions per minute. Lights continually signal rate progression, confirming compressions are within the 100- 120 cpm range recommended by the Industry

Guidelines. The yellow light will blink if the 120 cpm upper limit is exceeded, warning students to slow down. Instructors can monitor several students quickly and easily.

In addition to the visual feedback, students will also hear an audible beep as the chest is pushed to the appropriate depth of 2.0 inches allowing them to experience the true force needed to deliver real life chest compressions.

[Read More](#)



PRESTAN ADULT SERIES 2000 MANIKIN

SKU: PP-AM-2000-1-MS

The same high quality, realistic Professional Adult Manikin is now available with Advanced Feedback including Rate, Depth, Recoil, Ventilation, Hands-Off Time, and Chest Compression Fraction via a Bluetooth® enabled PRESTAN CPR Feedback App available on both Apple and Android devices.

PRESTAN PROFESSIONAL ADULT SERIES 2000 MANIKIN

The same high quality, realistic Professional Adult Manikin is now available with Advanced Feedback including Rate, Depth, Recoil, Ventilation, Hands-Off Time, and Chest Compression Fraction via a Bluetooth® enabled PRESTAN CPR Feedback App available on both Apple and Android devices.

Fully compliant with current industry guidelines, including the Heart and Stroke Foundation of Canada's February 2021 Use of Feedback Devices in Resuscitation Training Courses Policy and the most recent AHA Integrated Feedback Directive.

Available in medium or dark skin tones, with CPR Feedback only and in Single or 4-Pack configurations.

All PRESTAN CPR Training Manikins are Made in the USA.



ADDITIONAL DETAILS:

***REDESIGNED* PRESTAN CPR Feedback App**

Provides high quality, real time, advanced CPR feedback via a Bluetooth® connection with the ***NEW & IMPROVED*** PRESTAN Professional Adult Series 2000 Manikin.




Advanced CPR feedback (includes rate, depth, recoil, ventilation, hands-off time, and **chest compression fraction**) is monitored, easily displayed and available for reporting allowing students to self-monitor performance and gain confidence during CPR training.

The app allows instructors to provide comprehensive CPR training while monitoring CPR performance for up to six (6) students at a time and is equipped with a reporting feature that allows for CPR training results to be downloaded directly from the app.

The ***REDESIGNED*** PRESTAN CPR Feedback App provides an

- **improved user experience**
- **easier setup** - with no need to calibrate breaths
- more **intuitive CPR feedback** - with the addition of chest compression fraction
- **expanded list of available languages**

Languages currently supported by the PRESTAN CPR Feedback App include German, English (Intl), English (US), Spanish, French, Italian, Dutch, Russian, Korean, Arabic, Malay, Thai, Vietnamese, Belarusian, Croatian, Czech, Danish, Finnish, Greek, Hungarian, Norwegian, Polish, Portuguese, Romanian, Slovak, Slovene, Turkish, Ukrainian, and Chinese.



Download the ***REDESIGNED*** PRESTAN CPR Feedback App by Scanning the QR Code Below



The PRESTAN CPR Feedback App is for Instrumented Directive Feedback Devices (IDFD) that meet the most recent requirements of the American Heart Association and the Heart and Stroke Foundation of Canada's February 2021 *Use of Feedback Devices in Resuscitation Training Courses Policy*.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. FCC ID# QOQ-BGM220S2 / FCC ID# 2ASY3AK2G-PA

***NEW* PRESTAN Professional Adult Series 2000 Manikin CPR Rate Monitor LED Sequence**

	Connected Mode	Disconnected Mode
Sleep Mode/No Power		
Smart Device Connected Mode (Cycles in sequence left to right)		n/a
1 st Full Compression, Not Connected Mode (Wake mode/LED Test Mode)	n/a	
CPM* < 60		
CPM* 60 – 79		
CPM* 80 – 99		
CPM* 100 - 120		
CPM* > 120		
Compressions Stopped , Not Connected Mode		
Compressions Stopped , Not Connected Mode Low Battery Mode		
Compressions Stopped, Connected Mode		

*CPM = Compressions per minute
 = Flashing

In combination with PRESTAN’s traditional CPR feedback light system installed in the shoulder, the Series 2000 manikin offers a unique, dual-purpose training solution providing instructors with an option to train with either advanced or basic feedback.

The visual feedback gives students a real-life feel of delivering 100-120 compressions per minute. Lights continually signal rate progression, confirming compressions are within the 100- 120 cpm range recommended by the Industry Guidelines. The yellow light will blink if the 120 cpm upper limit is exceeded, warning students to slow down. Instructors can monitor several students quickly and easily.

In addition to the visual feedback, students will also hear an audible beep as the chest is pushed to the appropriate depth of 2.0 -2.4 inches allowing them to experience the true force needed to deliver real life chest compressions.

Series 2000 Add-On Kit

The PRESTAN Professional Adult Series 2000 Manikin Add-On Kit provides instructors with an option to upgrade their current PRESTAN Professional Adult and Jaw Thrust Manikins with the new Advanced Feedback parameters.

The PRESTAN Adult Series 2000 Add-On Kit is recommended for 2018 manikin models or newer.

Please Note: The newest version of the Series 2000 Add-On Kit contains a blue connector attached to the end of the white ventilation cord; however, the installation steps remain the same.

Series 2000 Lung Bags

The PRESTAN Professional Adult Series 2000 Lung Bags use the same easy to install PRESTAN design with the addition of a breath sensor. The Series 2000 lung bags are required for ventilation feedback via the PRESTAN CPR Feedback App.



[Read More](#)





PRESTAN ADULT MANIKINS


SKU: PP-AM-100M-MS

The PRESTAN Professional Training Manikin series is unlike any other on the market. Uniquely designed for realistic experiences, they are available individually or in sets of four for group training. CPR manikins are designed to match the real-life CPR experience.

Features:

- Fully compliant with current industry guidelines, including the most recent AHA Integrated Feedback Directive.
- Realistic looking and performing manikin is more acceptable to both instructors and students.
- Far easier and faster to carry and set up at training sites.
- Patented face/head tilt simulates how an actual victim's head would move if he required CPR and teaches students how to properly open the airway.
- Easily inserted face-shield/lung-bag provides visual chest rise to ensure students receive the most accurate and realistic training.
- The face shield is secured so it won't slip and there's no need for adhesives.
- Manikins are available with the CPR Feedback Monitor - Monitors may also be purchased separately to easily install into non-monitored Prestan Professional Manikins.
- Prestan CPR Manikins are not made with natural rubber latex.

Download PDF:



[Read More](#)



PRESTAN FEMALE MANIKIN

SKU: PP-AFM-100M-MS

The PRESTAN Professional Female Manikin is designed to transform CPR/AED training, improve CPR survival rates in women, encourage dialogue & breakdown barriers, and eliminate gender bias in CPR administration.

Using the same trusted design as the Professional Adult Manikin, the Female Manikin is realistic to the eye and the touch, unlike any other CPR Training Manikin on the market. Lightweight with fast and simple setup. Equipped with a revolutionary CPR Rate Monitor that allows for instant rate and depth feedback to both the instructor and student.

Fully compliant with current industry guidelines, including the most recent AHA Integrated Feedback Directive.

Available in medium or dark skin tones, with CPR Feedback only and in Single, 4-Pack or 4-Pack Diversity Kit configurations. Replacement skins are also available and can be used with the PRESTAN Professional Adult, Adult Series 2000, PRO+ Adult, and Adult Jaw Thrust Manikins.

All PRESTAN CPR Training Manikins are Made in the USA.



ADDITIONAL DETAILS:

Transform Training with Female CPR/AED Scenarios

When used in combination with the PRESTAN Female Accessory (to represent a bra), Manikin Shirt, and AED Trainer - PRESTAN offers the most realistic training experience to help build confidence when learning these critical lifesaving skills.



PRESTAN Professional Manikin CPR Rate Monitor


All PRESTAN Professional Manikins are available with a revolutionary CPR Rate Monitor that allows for instant feedback to both instructors and students. The visual feedback from the CPR Rate Monitor gives students a real life feel of delivering 100-120 compressions per minute. Lights continually signal rate progression, confirming compressions are within the 100- 120 cpm range recommended by the Industry Guidelines. The yellow light will blink if the 120 cpm upper limit is exceeded, warning students to slow down. Instructors can monitor several students quickly and easily.

In addition to the visual CPR Rate Monitor, students will also hear a clicker sound as the chest is pushed to the appropriate depth of 2.0 -2.4 inches allowing them to experience the true force needed to deliver real life chest compressions.

The CPR Rate Monitor may be purchased separately to upgrade any non-monitored PRESTAN Professional Manikin.

Sleep Mode/No Power	
1 st Full Compression (Wake mode/LED Test Mode)	
CPM* < 60	
CPM* 60 – 79	
CPM* 80 – 99	
CPM* 100 - 120	
CPM* > 120	
Compressions Stopped	

*CPM = Compressions per minute
 = Flashing



[Download PDF:](#)



Prestan Professional Manikins

ADULT MANIKIN

USE AND CARE INSTRUCTION SHEET

THANK YOU for buying the **Prestan Professional Adult Manikin!**

Please read through these instructions and keep them with your manikin(s). For information on Prestan's Authorized Distributors, help, guides, resources and more, visit the Prestan website at [www.prestan.asia](https://prestan.asia)



Prestan Professional Adult Manikin



Prestan Professional Adult Manikin

For future reference, make note below of the Prestan Authorized Distributor where you purchased your manikin. Call your Prestan Authorized Distributor for product questions, new device purchases, or warranty issues:

Prestan Authorized Distributor

Distributor Phone

Purchase Date



COPYRIGHT 2024 PRESTAN PRODUCTIONS LLC

[Read More](#)



PRESTAN PRO+ ADULT MANIKIN

SKU: PP-AMP-100M-MS

Using the same trusted design of the Professional Adult Manikin, PRESTAN is now introducing the PRO+ Adult Manikin featuring a new head design with functional oral and nasal airways along with an improved lung bag design for fast and easy installation.

With the introduction of the new PRO+ Adult Head design, PRESTAN is enhancing realism with a functional nasal airway and the ability to use a NARCAN(R) or NPA training device.

Fully compliant with current industry guidelines, including the most recent AHA Integrated Feedback Directive.

Available in medium or dark skin tones, with CPR Feedback only, in Single or 4-Pack configurations; the Professional Manikin Collection is available in Europe only.


The PRO+ Adult Head is also sold separately, allowing you to upgrade the Professional Adult Manikin.

All PRESTAN CPR Training Manikins are Made in the USA.

ADDITIONAL DETAILS:

Designed For Use With a NPA Device

When lubricated properly, a nasopharyngeal airway (NPA) can be inserted through the nose of the PRO+ Adult Head to simulate opening the nasal airway during training.



PRO+ Adult Lung Bags

The new PRO+ Adult Head design allows for a simplified lung bag installation; please refer to the PRO+ Adult Head Help Guide or Sales Video for more information on how to install the lung bag.



PRO+ Adult Face Shields

A new PRO+ Adult Face Shield has been designed to fully cover both oral and nasal airways; the current PRESTAN Professional Adult Face Shield can also be used.





PRESTAN Professional Manikin CPR Rate Monitor

All PRESTAN Professional Manikins are available with a revolutionary CPR Rate Monitor that allows for instant feedback to both instructors and students. The visual feedback from the CPR Rate Monitor gives students a real life feel of delivering 100-120 compressions per minute. Lights continually signal rate progression, confirming compressions are within the 100- 120 cpm range recommended by the Industry Guidelines. The yellow light will blink if the 120 cpm upper limit is exceeded, warning students to slow down. Instructors can monitor several students quickly and easily.


In addition to the visual CPR Rate Monitor, students will also hear a clicker sound as the chest is pushed to the appropriate depth of 2.0 -2.4 inches allowing them to experience the true force needed to deliver real life chest compressions.

The CPR Rate Monitor may be purchased separately to upgrade any non-monitored PRESTAN Professional Manikin.


Sleep Mode/No Power	
1 st Full Compression (Wake mode/LED Test Mode)	
CPM* < 60	
CPM* 60 – 79	
CPM* 80 – 99	
CPM* 100 - 120	
CPM* > 120	
Compressions Stopped	

*CPM = Compressions per minute

= Flashing



Download PDF:



[Read More](#)



PRESTAN INFANT ULTRALITE MANIKIN

SKU: #PP-IULM-100M

As with all PRESTAN manikins, the Infant Ultralite Manikin provides a realistic training experience with anatomically correct features, accurate head tilt and visible chest rise while providing intuitive feedback during CPR training.


Available in convenient and lightweight Single, 4-pack and 12-packs for efficient on-the-go training. This durable manikin offers an affordable solution for CPR training, with all the quality and realism you expect from PRESTAN.

Fully compliant with current industry guidelines, including the most recent AHA Integrated Feedback Directive.

Available in medium or dark skin tones in Single, 4-Pack or 12-Pack configurations; also available in 4-Pack or 12-Pack Diversity Kits.



Download PDF:



[Read More](#)



PRESTAN CHILD MANIKINS

SKU: PP-CM-100M-MS


The Prestan Professional Child Manikin is designed to look like a child verses an adult with less prominent body structures and softer, more childlike facial features. It is smaller and thinner than the Prestan Professional Adult Manikin. Yet, still realistic to the eye and the touch.

Features:

- Realistic looking and performing manikin is more acceptable to both instructors and students.
- Far easier and faster to carry and set up at training sites.
- Patented face/head tilt simulates how an actual victim's head would move if he required CPR and teaches students how to properly open the airway.
- Easily inserted face-shield/lung-bag provides visual chest rise to ensure students receive the most accurate and realistic training.
- The face shield is secured so it won't slip and there's no need for adhesives.
- Manikins are with the CPR Feedback Monitor
- Prestan Manikins are not made with natural rubber latex.



Download PDF:



[Read More](#)



PRESTAN INFANT MANIKIN

SKU: PP-IM-100M-MS


Life-like features and weight allows for realistic infant CPR training. Sculpted body contours and realistic skin feel help prepare students for real-life rescues. A head that simulates the way an actual infant's head would move. Face/head tilt enhances the CPR learning experience and teaches students how to open the airway during CPR. Available in single or 4 packs.

Features:

- Realistic looking and performing manikin is more acceptable to both instructors and students.
- Far easier and faster to carry and set up at training sites.
- Patented face/head tilt simulates how an actual victim's head would move if he required CPR and teaches students how to properly open the airway.
- Easily inserted face-shield/lung-bag provides visual chest rise to ensure students receive the most accurate and realistic training.
- The face shield is secured so it won't slip and there's no need for adhesives.
- Manikins are available with the CPR Feedback Monitor - Monitors may also be purchased separately to easily install into non-monitored Prestan Professional Manikins.
- Prestan Manikins are not made with natural rubber latex.



Download PDF:



[Read More](#)



PRESTAN ULTRALITE® MANIKINS


SKU: PP-ULM-100-MS

This is PRESTAN's most portable CPR manikin yet! It is amazingly easy to set up and sold in a convenient and lightweight package for efficient training-on-the-go. This durable vinyl and latex-free manikin offer an affordable tool for CPR training with all the quality and realism you expect from PRESTAN.

Features:

- Uniquely fast and easy to set-up.
- Stackable torso and head generate a compact package for transporting.
- Clamshell torso opening provides quick access for the easy-to-insert lung bag.
- Realistic skin resists dirt, grime, and grease.
- Adequate upper chest sizing facilitates AED pad placement and training.
- Visual anatomy landmarks reinforce proper hand placement for compressions.
- Head-tilt illustrates how to open the airway and allows each student to observe a true chest rise when the proper head-tilt chin-lift technique is used.
- Life tested to a half-million compressions.
- The manikin is easy to clean.

Download PDF:



[Read More](#)



PRESTAN COLLECTION CPR MANIKINS

SKU: PP-FM-300M-MS

The Prestan Professional Training Manikin series is unlike any other on the market. Uniquely designed for realistic experiences, they are available individually or in sets of four for group training. The manikins are designed to match the real-life CPR experience.



In CPR training, a precise feedback indicator is vital. Prestan's revolutionary CPR Rate Monitor provides instant feedback to both instructor and student on the correct compression depth. It enables each student to individually gauge their compression rate as well as allowing the instructor to monitor several students simultaneously. It is equipped with a unique Clicker Mechanism that gives an alert sound when the correct compression depth is reached. This helps ensure the students are using the correct force when administering CPR compressions.

Features:

- CPR feedback monitor that meets 2015 AHA CPR guidelines for Adult Manikin.
- Available in Adult, Child, and Infant models.
- Realistic looking and performing manikin is more acceptable to both instructors and students.
- Far easier and faster to carry and set up at training sites.
- Patented face/head tilt simulates how an actual victim's head would move if he required CPR and teaches students how to properly open the airway.
- Easily inserted face-shield/lung-bag provides visual chest rise to ensure students receive the most accurate and realistic training.
- The face shield is secured so it won't slip and there's no need for adhesives.
- Prestan Manikins are not made with natural rubber latex.



Download PDF:





Prestan Professional Manikins
INFANT MANIKIN
USE AND CARE INSTRUCTION SHEET

THANK YOU for buying the **Prestan Professional Infant Manikin!**

Please read through these instructions and keep them with your manikin(s). For information on Prestan's Authorized Distributors, help guides, resources and more, visit the Prestan website at www.prestan.asia/products/usa



Prestan Professional Infant Manikin



Prestan Professional Infant Manikin #400

For future reference, make note below of the Prestan Authorized Distributor where you purchased your manikin. Call your Prestan Authorized Distributor for product questions, reordering parts/accessories, or warranty issues.

Prestan Authorized Distributor

Distributor Phone

Purchase Date



[Read More](#)



<https://prestana.asia>

Page: 42

PRESTAN FAMILY PACK CPR MANIKINS

SKU: PP-FM-500M-MS





The Prestan Professional Training Manikin series is unlike any other on the market. Uniquely designed for realistic experiences, they are available individually or in sets of four for group training. The manikins are designed to match the real-life CPR experience.


In CPR training, a precise feedback indicator is vital. Prestan's revolutionary CPR Rate Monitor provides instant feedback to both instructor and student on the correct compression depth. It enables each student to individually gauge their compression rate as well as allowing the instructor to monitor several students simultaneously. It is equipped with a unique Clicker Mechanism that gives an alert sound when the correct compression depth is reached. This helps ensure the students are using the correct force when administering CPR compressions.

Features:

- CPR feedback monitor that meets the 2015 AHA CPR guidelines for Adult Manikin.
- Available in Adult, Child, and Infant models.
- Realistic looking and performing manikin is more acceptable to both instructors and students.
- Far easier and faster to carry and set up at training sites.
- Patented face/head tilt simulates how an actual victim's head would move if he required CPR and teaches students how to properly open the airway.
- Easily inserted face-shield/lung-bag provides visual chest rise to ensure students receive the most accurate and realistic training.
- The face shield is secured so it won't slip and there's no need for adhesives.
- Prestan Manikins are not made with natural rubber latex.

Download PDF:





[Read More](#)



PRESTAN ULTRALITE MANIKINS 4-PACK

SKU: PP-ULM-400-MS

The Prestan Professional Training Manikin series is unlike any other on the market. Uniquely designed for realistic experiences, they are available individually or in sets of four for group training. The manikins are designed to match the real-life CPR experience.


The Ultralite Manikin is amazingly easy to set up and is sold in a convenient and lightweight 4-pack for efficient training-on-the-go. This durable manikin offers an affordable tool for CPR training, with all the quality and realism you expect from Prestan.

Features:

- Stackable torsos and heads generate a compact package for transporting multiple units.
- The anatomically correct design includes the xiphoid process.
- Life tested to a half-million compressions.
- Visible chest rise with lung bag use.
- Compatible with AED Trainer pads.



Download PDF:



[Read More](#)



PRESTAN ADULT MANIKINS 4-PACK

SKU: PP-AM-400M-MS



The Prestan Professional Training Manikin series is unlike any other on the market. Uniquely designed for realistic experiences, they are available individually or in sets of four for group training. The manikins are designed to match the real-life CPR experience.

Realistic to the eye and the touch, the Prestan adult professional manikin is uniquely designed as a clamshell that accommodates an easy-to-insert face shield/lung bag.

Prestan's revolutionary new CPR Rate Monitor allows for instant feedback to both instructor and student regarding the rate of chest compression. This allows each student to gauge their rate of compressions on their own as well as allowing the instructor to monitor several students quickly and easily.


And with our patented face/head tilt, the Prestan Professional Adult Manikin simulates the way an actual victim's head would move if he required CPR.

In addition, the unique construction of our manikins incorporates a mechanism to help them use the correct force to compress the chest to the correct depth. It comes with 50 A/C face-shield lung bags and a nylon carrying case.

Features:

- Life-like features and weight allows for realistic adult/child CPR training
- Sculpted body contours and realistic skin feel help prepare students for real-life rescues
- A unique clamshell design that accommodates an easy-to-insert face shield lung bag
- CPR Rate Monitor allows for instant feedback to both instructor and student regarding the rate of chest compression
- Includes 4 Manikin with CPR Monitor, Face-Shield Lung Bags, and Nylon Carrying Case.

Download PDF:



[Read More](#)



PRESTAN INFANT MANIKINS 4-PACK

SKU: PP-IM-400M-MS

This AHA Compliant, 4-Pack of Prestan Infant Manikins-Medium Tone comes equipped with Compression Depth Clicker Feedback Devices and Compression Rate Monitor Feedback Devices.

This 4-Pack of Manikins will optimize your training program by ensuring the student's Compression Depth and Compression Rate are reached with an audible 'clicking' sound and LED indicators provided by 2 separate feedback devices. Additional Items Included: 50 Infant Lung Bags, an Instruction Sheet, and a nylon carrying case.

Achieve this with Prestan Professional's revolutionary new infant CPR manikin with Rate Monitor. The LED indicators allow for instant feedback to both instructors and students regarding their rate of chest compressions. Instructors can monitor several students quickly and easily, while students themselves can visually gauge their compression rate reinforcing more practice time using proper technique.


The depth of compression is a crucial skill to learn as well. Class participants learn proper compression depth with the manikin's audible 'clicking' sound. When compressions are done to the ECC recommended depth, the manikin's 'click' lets students know they are on target. Compression rate and depth – two vitally important components to learning for successful CPR. With Prestan's visual and audible reinforcement of good CPR techniques, students will confidently embrace their role in the cardiac chain of survival.

Features:

- LED indicators on the infant's diaper allow for instant feedback to both instructor and student regarding the accuracy of chest compressions.
- The realistically sized 100% Latex Free infant has a pinch-able nose. "Flick their feet" reminds students to check for the level of consciousness – easily taught since the manikin has fully formed legs with feet. Moveable arms even allow for brachial pulse checks.
- Practice conscious choking back blows on the infants' padded upper back portion.
- The easy to insert lung-bag speeds manikin setup.
- The face-shield is secured so it won't slip and there's no need for adhesives.
- Realistic skin resists dirt, grime, grease, and even adhesive build up from AED electrodes!
- Visual anatomy 'landmarks' such as clavicles, an accurate nipple line, sternum, and xiphoid process facilitates the practice of choking skills as well as reinforcing proper finger placement for compressions.
- A tiltable head illustrates how to open the airway and allows each student to observe a true chest rise when a proper head-tilt chin-lift technique is used.
- The infant's head structure is capable of demonstrating the jaw-thrust maneuver.
- Audible clicker mechanism increases the student's confidence in administering adequate depth of compressions.
- Manikins are capable of being cleaned and decontaminated per CDC guidelines.



Download PDF:



[Read More](#)



PRESTAN INSTRUCTOR KIT

SKU: PP-INSTRKIT-1-MSDS

The PRESTAN Instructor Kit combines our most popular CPR training manikins (two Professional Adults w/ Female Accessory & two Professional Infants) along with two AED UltraTrainers® and a collection of the most commonly used CPR Training Supplies in one single collection, packaged in a heavy-duty canvas bag with wheels & telescoping handle for enhanced portability.

Fully compliant with current industry guidelines, including the most recent AHA Integrated Feedback Directive.

Available as a Diversity Kit (two medium skin tone and two dark skin tone manikins) and with CPR Feedback only.

All PRESTAN CPR Training Manikins are Made in the USA.

ADDITIONAL DETAILS:

PRESTAN Instructor Kit Detailed Product Contents

- (2) Professional Adult Manikin w/ CPR feedback (1 Medium Skin & 1 Dark Skin)
- (2) Professional Infant Manikin w/ CPR feedback (1 Medium Skin & 1 Dark Skin)
- (2) AED UltraTrainer® (English/Spanish)
- (2) Female Accessory (1 Medium Skin & 1 Dark Skin)
- (2) Adult CPR Manikin Training Shirt
- (2) Infant CPR Manikin Training Shirt
- (50) Adult Face-Shield/Lung-Bag
- (50) Infant Face-Shield/Lung-Bag
- (2) Adult Bag Valve Mask
- (2) Infant Bag Valve Mask
- (4) PRESTAN Kneeling Pad
- (1) Carry Bag w/ Roller Wheels & Telescopic Handle


PRESTAN Professional Manikin CPR Rate Monitor

All PRESTAN Professional Manikins are available with a revolutionary CPR Rate Monitor that allows for instant feedback to both instructors and students. The visual feedback from the CPR Rate Monitor gives students a real life feel of delivering 100-120 compressions per minute. Lights continually signal rate progression, confirming compressions are within the 100- 120 cpm range recommended by the Industry Guidelines. The yellow light will blink if the 120 cpm upper limit is exceeded, warning students to slow down. Instructors can monitor several students quickly and easily.


In addition to the visual CPR Rate Monitor, students will also hear a clicker sound as the chest is pushed to the appropriate depth of 2.0 -2.4 inches allowing them to experience the true force needed to deliver real life chest compressions.

The CPR Rate Monitor may be purchased separately to upgrade any non-monitored PRESTAN Professional Manikin.



Sleep Mode/No Power	
1 st Full Compression (Wake mode/LED Test Mode)	
CPM* < 60	
CPM* 60 – 79	
CPM* 80 – 99	
CPM* 100 - 120	
CPM* > 120	
Compressions Stopped	

*CPM = Compressions per minute
 = Flashing



[Download PDF:](#)



Prestan Professional Manikins ADULT MANIKIN

USE AND CARE INSTRUCTION SHEET

THANK YOU for buying the **Prestan Professional Adult Manikin!**
Please read through these instructions and keep them with your manikin(s). For information on Prestan's Authorized Distributors, help guides, resources and more, visit the Prestan website at www.prestanproducts.com



Prestan Professional Adult Manikin Prestan Professional Adult Manikin 4-Pk

For future reference, make note below of the Prestan Authorized Distributor where you purchased your manikin. Call your Prestan Authorized Distributor for product questions, ordering supplies/accessories, or warranty issues.

Prestan Authorized Distributor	Distributor Phone	Purchase Date

© 2016 PRESTAN PRODUCTS LLC





Prestan Professional Manikins INFANT MANIKIN

USE AND CARE INSTRUCTION SHEET

THANK YOU for buying the **Prestan Professional Infant Manikin!**
Please read through these instructions and keep them with your manikin(s). For information on Prestan's Authorized Distributors, help guides, resources and more, visit the Prestan website at www.prestanproducts.com



Prestan Professional Infant Manikin Prestan Professional Infant Manikin 4-Pk

For future reference, make note below of the Prestan Authorized Distributor where you purchased your manikin. Call your Prestan Authorized Distributor for product questions, ordering supplies/accessories, or warranty issues.

Prestan Authorized Distributor	Distributor Phone	Purchase Date

© 2016 PRESTAN PRODUCTS LLC



Prestan® AED UltraTrainer™





PRESTAN®



[Read More](#)



PRESTAN PROFESSIONAL ADULT 50PCS FACE-SHIELD / LUNG BAG


SKU: PP-ALB-50

Manikin Replacement Lung Bags

Features:

- Prestan Manikin Replacement Lung Bags.
- Clear plastic inflatable lungs.
- Easy to insert.
- Disposable.

[Read More](#)





PRESTAN PROFESSIONAL CHILD 50PCS FACE-SHIELD

SKU: PP-CFS-50

In classes where there is more than one student per manikin, there is a need for a balance between protecting the student and giving the student a realistic experience.

This is a 50 pack of Filtered Face Shields for Prestan Professional CPR Manikins. If manikins are shared during your training class, you'll want to use a filtered face shield to minimize the risk of cross-contamination between students. The Prestan face shield makes it easy and since it's secured to the manikin, it won't slip. No need for adhesives! Your students will spend less time distracted by an unruly barrier and more time practicing their skills.

Every CPR/AED instructor wants the best outcome when and if the student has a real-life experience giving CPR. This face shield, placed over the manikin's face, allows each student to practice CPR, be protected, and have the benefit of a realistic experience of a full chest rise.

50 pack face shields for Prestan Child Manikins.

Features:

- Professional Child Face Shields by Prestan Products (50 Pack)
- Not made with natural rubber latex

[Read More](#)



PRESTAN PROFESSIONAL CHILD 50PCS FACE-SHIELD / LUNG BAG

SKU: PP-CLB-50

Replacement for Manikin Face Shield Lung Bags.

Features:

- Prestan Manikin Replacement Face Shield Lung Bags
- Clear plastic inflatable lungs
- Easy to insert
- Disposable

[Read More](#)



PRESTAN PROFESSIONAL INFANT 50PCS FACE-SHIELD / LUNG BAG

SKU: PP-ILB-50

Replacement for Manikin Face Shield Lung Bags

Features:

- Prestan Manikin Replacement Face Shield Lung Bags
- Clear plastic inflatable lungs
- Easy to insert
- Disposable

[Read More](#)



PRESTAN PROFESSIONAL ULTRALITE 50PCS FACE-SHIELD / LUNG BAG

SKU: PP-ULB-50

The face shield lung bag delivers a realistic experience of a full chest rise during CPR training. These lung bags fit the Prestan Ultralite series CPR manikins only.

This pack includes 50 lung bags. Insertion of the one-piece lung bag design is a breeze! Fold it up, slide it in and you are ready for training.

Features:

- Easy to insert
- Work with Prestan Ultralite CPR Manikins
- No tools necessary to insert

[Read More](#)



PRESTAN FEMALE ACCESSORY

SKU: PP-ULB-50-1




The PRESTAN Female Accessory™ (*patent-pending*) has been designed to build confidence in CPR training, improve CPR survival rates in women, encourage dialogue & breakdown barriers, and eliminate gender bias in CPR administration.

Securely fits both the PRESTAN Professional Adult/Series 2000 & Ultralite Manikins. Pair with the PRESTAN CPR Training Shirts for a full real-life training experience.

Available in packs of 4 and in medium or dark skin tones.

[Read More](#)



PRESTAN MASKS ADAPTERS

SKU: 10076-PPA

Use this affordable, disposable adapter in place of expensive mask valves for cost-effective training. It features an integral air filter that obstructs the flow of contaminants. Fits all standard rescue masks (15mm I.D. or 22mm O.D.).

***For training purposes only.**

Features:

- It features an integral air filter that obstructs the flow of contaminants.
- Fits all standard rescue masks
- Rescue mask not included

[Read More](#)





PRESTAN AED ULTRA TRAINER

SKU: PP-AEDUT-101

The AED for training purposes. PRESTAN is the First and Only Universal AED Trainer with **2015 Guidelines Compliant and Upgradable**.

The PRESTAN AED UltraTrainer is equipped with a dual graphics pad and Prestan's unique Pad Sensing System, which automatically detects when each pad is placed on the manikin. This AED trainer can accurately simulate a live AED, allowing for a more realistic training class.

The pads are durable, lasting up to 100 applications, and do not leave behind that annoying sticky residue. Replacement is easy and affordable.

The PRESTAN AED UltraTrainer is the first and only universal trainer to have a child button that best simulates modern live AEDs, giving your students a realistic training experience.

The easy to use audio menu allows you to customize your training with options of automatic or semi-automatic shock, 5 training scenarios, and the ability to turn on or off the compression metronome and ventilation command.

Features:

- PRESTAN is the First and Only Universal AED Trainer with **2015 Guidelines Compliant and Upgradable**
- **Play/Pause Button**
- **5 Training Scenarios**
- Dual Graphic (Adult/Child) Training Pads
- Pad Sensing System
- Preconnected Pads
- Child Training Button
- Child AED specific voice prompts
- Automatic or Semi-Automatic shock simulation
- CPR Coaching, "Give Breaths" voice prompt
- Turn on or off compression metronome
- Turn on or off "Give Breaths" prompts
- Low Battery warning

Download PDF:



[Read More](#)



PRESTAN PROFESSIONAL AED TRAINER PLUS

SKU: PP-AEDT2-101

QUALITY / DURABILITY / MOST SIMILAR TO A LIVE AED

- Leading-edge AED Trainer that most closely resembles brand AED Trainers in size and weight for a realistic training experience.
- It features customizable options to match training environments and includes a child button.
- Clear confident voice prompts. Also, CPR prompts to help the student keep pace with chest compressions.
- Easy to replace modules containing current CPR guidelines and 2 languages.
- A remote is optional and available for instructors who prefer this feature. A single remote will control several AED Trainers within its range.



Change the Language with a simple press of a button



Scenario button can quickly change the shock sequence (no need for manual)



Three volume options: loud, medium, soft



Play and Pause button for those teachable moments

Features:

- **Compliant with current 2015 AHA Guidelines**
- Fast and Easy as well as Update
- Adult/Child selection buttons
- Training pads with dual adult/child graphics
- Brightly illuminated shock button that flashes when the shock is advised
- 5 preprogrammed training scenarios with Clear and calm voice prompts
- Chest compression rate metronome with audible count reminders
- Automatic or semi-automatic shock delivery
- Programming option to either include or not include time for breaths after the compression cycle
- Available remote control (**sold separately**)
- Long battery life with automatic shut-off
- Replaceable language module
- Sturdy, flip-up lid

Download PDF:



[Read More](#)





PRESTAN PROFESSIONAL AED TRAINER REMOTE

SKU: PP-AEDT-100-R

To use with the Prestan Professional AED trainer. Remote for Prestan Professional AED Trainer by using 2 AA batteries. Remote control lets the Instructor select a scenario, pause to give instruction, and then resume the scenario in progress.

Features:

- Set-up buttons
- Unit control buttons:
- Simulate loose pad
- Initiate "Press Deeper" during CPR
- Pause
- Play
- Send scenario

[Read More](#)



PRESTAN PROFESSIONAL AED TRAINER PLUS PADS

SKU: PP-APAD2-4

PRESTAN Professional's high-performance training pads feature clear anatomical placement graphics and a long-life adhesive.

The pads are made of the same foam material used in live AED pads and conform well to all CPR manikins on the market today. The pad's silicone-based adhesive, coupled with their foam structure, ensure easy application and removal from your manikin.


The manufacturer states that they are good for at least 25 sticks and release AED training skill sessions, but we've used them for many more than that! The AED training pads ship in a zippered plastic bag, but once connected, stay clean and protected in the unit's training pad case.


***Pads are for training purposes only and cannot be used to deliver a defibrillation shock.**

Features:

- Fits Prestan AED Trainers
- Includes pad sensor
- Easy and clean removal from manikins
- Works with any brand manikin
- Latex-free silicone adhesive

[Read More](#)





PRESTAN AED ULTRATRAINER PADS (4 PAIRS IN A PACK)

SKU: PP-UTPAD-4

Prestan AED UltraTrainer™ Adult/Child Replacement Training Pad Set, 4-Pack (8 pads total).

Prestan AED UltraTrainer pads are durable and long-lasting and feature an exclusive pad sensing system that allows the trainer to determine the pads are attached to the manikin. The silicone-based adhesive provides easy application and removal, helping prolong the life of the manikin skin.

Features:

- Dual Graphic (Adult/Child) Training Pads
- Pad Sensing System
- Preconnected Pads

[Read More](#)



INDEX

P

PRESTAN Adult Manikins 29
PRESTAN Adult Manikins 4-Pack 46
PRESTAN Adult Series 2000 Manikin 25
PRESTAN AED Ultra Trainer 60
PRESTAN AED UltraTrainer Pads (4 Pairs In A Pack) 10, 66
PRESTAN Child Manikins 38
PRESTAN Child Series 2000 Manikin 21
PRESTAN Collection CPR Manikins 41
PRESTAN Family Pack CPR Manikins 43
PRESTAN Female Accessory 7, 58
PRESTAN Female Manikin 30
PRESTAN Female Series 2000 Manikin 12
PRESTAN Infant Manikin 39
PRESTAN Infant Manikins 4-Pack 48
PRESTAN Infant Series 2000 Manikin 17
PRESTAN Infant Ultralite Manikin 37
PRESTAN Instructor Kit 50
PRESTAN Masks Adapters 8, 59
PRESTAN PRO+ Adult Manikin 33
PRESTAN Professional Adult 50pcs Face-Shield / Lung Bag 6, 53
PRESTAN Professional AED Trainer PLUS 62
PRESTAN Professional AED Trainer PLUS Pads 11, 65
PRESTAN Professional AED Trainer Remote 9, 64
PRESTAN Professional Child 50pcs Face-Shield 5, 54
PRESTAN Professional Child 50pcs Face-Shield / Lung Bag 4, 55
PRESTAN Professional Infant 50pcs Face-Shield / Lung Bag 3, 56
PRESTAN Professional Ultralite 50pcs Face-Shield / Lung Bag 2, 57
PRESTAN Ultralite Manikins 4-Pack 45
PRESTAN Ultralite® Manikins 40